



Ocular Four Mirror Mini Gonio Lenses									
CE	Product Code	Image Mag	Laser Spot Mag	Contact OD	Lens Height	Ring OD	Static Gonio FOV		
	ARGON/DIODE LASER							<i>Journal Reference: Optometric Management, Vol. 35, No. 6, June 2000</i>	
	O4GFA	.94x	1.06x	15mm	23.5mm	23.5mm	120°		
	O4GFA-LR	.94x	1.06x	15mm	26.9mm	32.3mm	120°		
	DIAGNOSTIC								
	O4GF	.94x	na	15mm	22.5mm	23.5mm	120°		
O4GF-LR	.94x	na	15mm	25.8mm	32.3mm	120°			

Lens Design

- For anterior chamber observation and photocoagulation procedures.
- The Four Mirror Mini Gonio Lens was developed in conjunction with Asian doctors to increase the ease of use and examination of Asian eyes and small palpebral fissures.
- The small diameter endpoint allows the lens to be tilted slightly in either direction for optimum viewing and makes it ideal for use on children or adults with small palpebral fissures.
- Its four mirrors are inclined at 62° and are positioned 90° apart to allow complete observation of angle with little rotation.
- Field of View through central window = 33°
- A broadband anti-reflective coating is bonded to the O4GFA and O4GFA-LR lenses to minimize reflections and maximize light transmission during Argon/Diode Laser Treatment.
- Available with a standard or large holding ring.

Technique

- After the lens is placed on the anesthetized eye, indirect observation is used.
 - With the mirror placed at 12:00 using a narrow slit beam at approximately 10°, a section of the angle can be observed at the 6:00 area.
 - To observe the 3:00 and 9:00 areas, the slit lamp should be rotated in a horizontal position.
- Four mirrors allow the lens to only need a slight rotation to view the entire anterior chamber angle.
- The mirrors may be set at either of two positions (the square position or the diamond position).
 - The latter position, with the mirrors at 1:30, 4:30, 7:30 and 10:30 meridians, permits the slit lamp beam to be readily used in all four quadrants of the angle.
- Examination of a narrow angle can be facilitated in two ways.
 - The lens can be shifted or rocked slightly on the corneal surface in any direction for a millimeter or so.
 - This often brings a hidden angle recess or structure into view.
 - This can be further aided by shifting the position of the fixation light in the required direction.

Caution	When using the lens for photocoagulation, use extreme care to keep the laser away from the edges. If the beam strikes the area around the mirror, it may be absorbed and burn the area. Mirrors damaged in this manner cannot be repaired.
----------------	--

Cleaning

- Rinse: Immediately upon removal from patient's eye, thoroughly rinse in cool or tepid water.
- Wash: Place a few drops of mild soap on a moistened cotton ball. Gently clean with a circular motion.
- Rinse: Thoroughly rinse in cool or tepid water, then dry carefully with a *non-linting* tissue.
- Then: Proceed with either disinfection or sterilization instructions.

Disinfection										
Soak In:	<table border="1"> <thead> <tr> <th>GLUTARALDEHYDE</th> <th>OR</th> <th>BLEACH</th> </tr> </thead> <tbody> <tr> <td>2% or 3.4% aqueous solution</td> <td rowspan="3"></td> <td>10% solution mixed at: 1 part bleach to 9 parts cool tepid water</td> </tr> <tr> <td>Temperature per manufacturer instructions</td> <td rowspan="2">Recommended exposure time = 10 minutes</td> </tr> <tr> <td>Minimum exposure time = 20 minutes</td> </tr> </tbody> </table>	GLUTARALDEHYDE	OR	BLEACH	2% or 3.4% aqueous solution		10% solution mixed at: 1 part bleach to 9 parts cool tepid water	Temperature per manufacturer instructions	Recommended exposure time = 10 minutes	Minimum exposure time = 20 minutes
	GLUTARALDEHYDE	OR	BLEACH							
	2% or 3.4% aqueous solution		10% solution mixed at: 1 part bleach to 9 parts cool tepid water							
	Temperature per manufacturer instructions		Recommended exposure time = 10 minutes							
Minimum exposure time = 20 minutes										
CAUTION <i>To avoid damage to the lens, do not exceed recommended exposure time.</i>										
Then:	Rinse lens <i>thoroughly</i> to remove disinfection solution. 3 cycles of 1 minute, with cool or tepid water is recommended. Dry carefully and place in a dry storage case.									
NOTE	This lens is known to be compatible with: Asepti-Wipe, Cavi-cide, Cidex, Cidex OPA, DisCide Wipe, Enviro-cide, H ₂ O ₂ - 3%, and Opti-Cide									
CAUTION	<i>If used on an ulcerated cornea, lens must be STERILIZED before next procedure.</i>									

Sterilization						
AUTOCLAVE	STERRAD	STERIS SYSTEM 1	ETO	ETO Parameters		
No	No	YES	YES	Minimum Time	Temperature	Aeration Time
		Per manufacturer instructions	See Right	1 hour	130°F (54°C)	12 hours
WARNING	<i>Never Steam Autoclave or Boil listed lenses. Never soak in Alcohol, Acetone or Other Solvents.</i>					

For information on compatibility with alternative product care methods, contact Customer Service.

