




Ocular Fundus Lenses									
CE	Product Code	Style ⁽¹⁾	Image Mag	Laser Spot Mag	Contact OD	Lens Height	Static FOV		
	ARGON/DIODE								
	OGFA		.93x	1.08x	15mm	16.5mm	36°		
	OGFA-2	NMR	.93x	1.02x	15mm	16.5mm	36°		
	DIAGNOSTIC								
	OGF		.93x	na	15mm	16.5mm	36°		
	OGF-2	NMR	.97x	na	15mm	16.5mm	35°		

⁽¹⁾NMR "No Methylcellulose Required" feature was developed with Frederick M. Kapetansky, M.D., Columbus, Ohio

Lens Design

- The Fundus Diagnostic and Laser Lenses provide clear visualization of the posterior pole.
- The surface in contact with the cornea is contoured to match the corneal curvature and eliminate air bubbles.
- The lenses will stay in place without holding or may be manipulated by the cone.
- A broad-band, anti-reflective coating for Argon and Diode lasers is standard.
- An unusually flat cornea (K=38.00) may require use of a drop of methylcellulose or Celluvisc between the cornea and the lens on the NMR styles.

Cleaning

- Rinse: Immediately upon removal from patient's eye, thoroughly rinse in cool or tepid water.
Wash: Place a few drops of mild soap on a moistened cotton ball. Gently clean with a circular motion.
Rinse: Thoroughly rinse in cool or tepid water, then dry carefully with a *non-linting* tissue.
Then: Proceed with either disinfection or sterilization instructions.

Disinfection

Soak In:	GLUTARALDEHYDE		OR	BLEACH	
	2% or 3.4% aqueous solution			10% solution mixed at: 1 part bleach to 9 parts cool tepid water	
	Temperature per manufacturer instructions			Recommended exposure time = 10 minutes	
	Minimum exposure time = 20 minutes				
Caution To avoid damage to the lens, do not exceed recommended exposure time.					
Then:	Rinse lens <i>thoroughly</i> to remove disinfection solution. 3 cycles of 1 minute, with cool or tepid water is recommended. Dry carefully and place in a dry storage case.				
NOTE	This lens is known to be compatible with: Asepti-Wipe, Cavi-cide, Cidex, Cidex OPA, DisCide Wipe, Enviro-cide, H ₂ O ₂ - 3%, and Opti-Cide				
Caution	If used on an ulcerated cornea, lens must be <i>STERILIZED</i> before next procedure.				

Sterilization

AUTOCLAVE	STERRAD	STERIS SYSTEM 1	ETO	ETO Parameters		
No	No	YES	YES	Minimum Time	Temperature	Aeration Time
		Per manufacturer instructions	See Right	1 hour	130°F (54°C)	12 hours

WARNING Never Steam Autoclave or Boil listed lenses.
Never soak in Alcohol, Acetone or Other Solvents.

For information on compatibility with alternative product care methods, contact Customer Service.

