


## Ocular Hoskins-Barkan Goniotomy Lenses

	<b>Product Code</b>	<b>Style</b>	<b>Size</b>	<i>Designed with: H. Dunbar Hoskins, M.D San Francisco, CA</i>
	<b>OHBG-1</b>	Infant	10mm	
	<b>OHBG-2</b>	Premature Infant	9mm	
	<b>OHBG-3</b>	Adult	11.5 mm	
<b>CE</b>				

### Lens Design

- The Hoskins-Barkan Goniotomy Lens was primarily designed for transverse goniotomy surgery with the operating microscope, but can also be used as a diagnostic lens.
- The infant lens is oval and conical in shape, with a 10mm diameter magnified view of the anterior chamber and anterior chamber angle.
- The premature infant lens is the same in shape and design except the dimensions are 1mm smaller for premature infant surgery.
- There is also an adult size at 11.5 mm diameter.

### Cleaning

Rinse: Immediately upon removal from patient's eye, thoroughly rinse in cool or tepid water.  
 Wash: Place a few drops of mild soap on a moistened cotton ball. Gently clean with a circular motion.  
 Rinse: Thoroughly rinse in cool or tepid water, then dry carefully with a *non-linting* tissue.  
 Then: Proceed with either disinfection or sterilization instructions.

### Disinfection

Soak In:	<b>GLUTARALDEHYDE</b>	<b>OR</b>	<b>BLEACH</b>
	2% or 3.4% aqueous solution		10% solution mixed at: 1 part bleach to 9 parts cool tepid water
	Temperature per manufacturer instructions		Recommended exposure time = 10 minutes
	Minimum exposure time = 20 minutes		
<b>Caution</b> <i>To avoid damage to the lens, do not exceed recommended exposure time.</i>			
Then:	Rinse lens <i>thoroughly</i> to remove disinfection solution. 3 cycles of 1 minute, with cool or tepid water is recommended. Dry carefully and place in a dry storage case.		
<b>NOTE</b>	This lens is known to be compatible with: Asepti-Wipe, Cavi-cide, Cidex, Cidex OPA, DisCide Wipe, Enviro-cide, H <sub>2</sub> O <sub>2</sub> - 3%, and Opti-Cide		

### Sterilization

AUTOCLAVE	STERRAD	STERIS SYSTEM 1	ETO	ETO Parameters		
No	No	Yes	Yes	Minimum Time	Temperature	Aeration Time
		Per manufacturer instructions	See Right	1 hour	130°F (54°C)	12 hours
<b>WARNING</b>		<b>Never Steam Autoclave or Boil listed lenses. Never soak in Alcohol, Acetone or Other Solvents.</b>				

For information on compatibility with alternative product care methods, contact Customer Service.

