




Ocular 132D Indirect Vitrectomy Lens

	Product Code	Image Mag	Static FOV	Dynamic FOV	
	OIV-132 	.45x	100°	135°	

Design

Designed to be used on the Ocular Landers Wide Angle Surgical Viewing System (OSVS) in conjunction with an Inverter Vitrectomy System.

Attaches to the OSVS with the Ocular 132D Indirect Vitrectomy Lens Holder (OIV-H132)

Non-contact design allows patient's eye to be rotated freely to view the peripheral retina and vitreous

Steam sterilizable; sterilizing case included

Designed to allow a clear view in the fluid or air filled eye

Cleaning

Rinse: Immediately upon removal from patient's eye, thoroughly rinse in cool or tepid water.

Wash: Place a few drops of mild soap on a moistened cotton ball. Gently clean with a circular motion.

Rinse: Thoroughly rinse in cool or tepid water, then dry carefully with a *non-linting* tissue.

Then: Proceed with either disinfection or sterilization instructions.

Disinfection

Soak In:	GLUTARALDEHYDE	OR	BLEACH
	2% or 3.4% aqueous solution		10% solution mixed at: 1 part bleach to 9 parts cool tepid water
	Temperature per manufacturer instructions		Recommended exposure time = 10 minutes
	Minimum exposure time = 20 minutes		
Caution <i>To avoid damage to the lens, do not exceed recommended exposure time.</i>			
Then:	Rinse lens <i>thoroughly</i> to remove disinfection solution. 3 cycles of 1 minute, with cool or tepid water is recommended. Dry carefully and place in a dry storage case.		
NOTE	This lens is known to be compatible with: Asepti-Wipe, Cavi-cide, Cidex, Cidex OPA, DisCide Wipe, Enviro-cide, and Opti-Cide.		

Sterilization - Autoclave

Prep: Place cleaned lens on its side in sterilization case.

Process:	Standard Cycle (wrapped)				
	Temperature	Time	OR	Temperature	Time
	270°F (134°C)	15 minutes minimum		250°F (121°C)	30 minutes minimum
Caution	<i>Use only distilled water in the steam sterilizer.</i> If not distilled, mineral deposits from hard water (steam) will leave a cloudy film on the lens. The deposit can only be removed by regrinding and re-polishing the lens and repair costs approximate that of a new lens.				

Store: Place in a biological barrier peel pouch to ensure sterility after the process.

For Immediate use only	Flash autoclave (unwrapped) at a minimum of 270°F (134°C) for a minimum of 10 minutes.
-------------------------------	--

ETO		
Minimum Time	Temperature	Aeration Time
1 hour	130°F (54°C)	12 hours

Steris System 1
No

Sterrad
No

For information on compatibility with alternative product care methods, contact Customer Service.

