OCULAR INSTRUMENTS

2255 116th Ave NE, Bellevue, Washington 98004-3039 USA

T: 425-455-5200 or 800-888-6616 F: 425-462-6669

E: ocular@ocular-instruments.com I: www.ocular-instruments.com

Ocular Tonometers

OBVT
Barraquer Varley 90mm Hg

OBT-TC-10-15
Barraquer Terry 10-15mm Hg

OBBT
Baron Barraquer 65/90mm Hg

OBT-TC-15-21
Barraquer Terry 15-21mm Hg

OBT65/90
Barraquer 65/90mm Hg

OKBT-20/30
Kasaby Barraquer 20/30mm Hg

OBT-65
Barraquer 65mm Hg

OGBT-30/50
Griffin Barraquer 30/50mm Hg

Ocular Tonometers (English)

The engraved ring reticles on the endpoint indicates a predetermined intraocular pressure: **OBVT**: 90mm Hg. **OBBT**: 65mm Hg or 90mm Hg, the smaller ring is 90mm Hg. **OBT65/90**: 65mm Hg or 90mm Hg. The smaller ring is 90mm Hg. **OBT-65**: 65mm Hg. **OBT-TC-10-15**: 10mm Hg or 15mm Hg. The smaller ring is 15mm Hg. **OBT-TC-15-21**: 15mm Hg or 21mm Hg. The smaller ring is 21mm Hg. **OKBT-20/30**: 20mm Hg or 30mm Hg. The smaller ring is 30mm Hg. **OGBT-30/50**: 30mm Hg or 50mm Hg. The smaller ring is 30mm Hg. *Also Available for all tonometers: Silicone Replacement Ring OBT-O (5 pack)*

Lens Design: The Tonometers measure the intraocular pressure. They are based on Maklakov’s principle of applanation tonometry and consists of a conical shaft of high-density transparent plastic, a silicone retaining ring and a slip ring handle. **Technique**: The eye is pressurized until the intraocular pressure appears to be at the desired level. Any excessive fluid on the cornea is dried. The tonometer is gently lowered onto the cornea until the applanation circle can be seen through the tonometer. If the applanation circle is smaller than the inside of the engraved ring, pressure in the eye is lower than the calibrated value.

Cleaning: Rinse immediately upon removal from patient’s eye, thoroughly rinse in cool or tepid water. Wash by placing few drops of mild soap on a moistened cotton ball. Gently clean with a circular motion. Rinse thoroughly in cool or tepid water, and then dry carefully with a *non-linting* tissue. Then proceed with either disinfection or sterilization instructions.

Disinfecting: Tonometers can be soaked in GLUTARALDEHYDE (2% or 3.4% aqueous solution, temperature per manufacturer instructions with minimum exposure of 20 mins) OR BLEACH (10% solution mixed at 1 part bleach to 9 parts cool or tepid water, recommended exposure time is 10 mins). **Then** Rinse lens *thoroughly* to remove disinfection solution.3 cycles of 1 minute, with cool or tepid water is recommended. Dry carefully and place in a dry storage case. *Note*: This lens is known to be compatible with Asepti-Wipe, Cavi-cide, Cidex OPA,DisCide Wipe, Enviro-cide H₂O₂-3%, and Opti-Cide.

Sterilizing: Suitable for autoclaving. Only Flash for 4 minutes. Not compatible with Sterrad. Compatible with Steris System 1, per Manufacturer’s instructions. Compatible EO (EO parameters: minimum time 1 hour at temperature 130°F or 54° C and aeration time of 12 hours. **WARNING – Never soak in acetone, alcohol or other solvents. Autoclave – Flash Only**: Disassemble and thoroughly wash the tonometer so that it is free of mucous, sebaceous deposits or other debris. Place all three parts in a tray taking care to protect the tonometer from damage by other instruments. **Flash Autoclave Only** (unwrapped) for four (4) minutes at 270°F or 132°C. No dry time. **WARNING: REMOVE PROMPTLY, longer exposures will damage lens. The intense heat for an extended time will cause the plastic to cloud. Reassemble before use. In the absence of the ring, a false reading will occur.**

Tonometers have a lifetime of 5 years. After a period of 2 years of purchase, check for the following: Any visual damage, easy gliding and turning without any resistance, no complete rip of the white ‘O’ type joint ring, scratches on applanation (contact surface), complete visibility of engraved white ring on the applanation (contact) surface. Contact customer service for any such concern. Also for information on compatibility with alternative sterilization methods, contact Customer Service.

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