




## Ocular Kaufman Single Mirror Laser Lens

	Product Code	Contact OD	Lens Height	
	OKSMA 	13mm	19.3mm	

### Lens Design

- A single mirror, set at 64°, is designed for visualization and laser procedures of the anterior chamber in all species of monkey.
- Ocular's Laserlight® high efficiency, broad band, anti-reflective coating provides optimal image contrast, minimizes bothersome reflections and maximizes visible and diode laser transmission.

### Cleaning

Rinse: Immediately upon removal from patient's eye, thoroughly rinse in cool or tepid water.  
 Wash: Place a few drops of mild soap on a moistened cotton ball. Gently clean with a circular motion.  
 Rinse: Thoroughly rinse in cool or tepid water, then dry carefully with a *non-linting* tissue.  
 Then: Proceed with either disinfection or sterilization instructions.

### Disinfection

Soak In:	<b>GLUTARALDEHYDE</b>		<b>OR</b>	<b>BLEACH</b>		
	2% or 3.4% aqueous solution			10% solution mixed at: 1 part bleach to 9 parts cool tepid water		
	Temperature per manufacturer instructions			Recommended exposure time = 10 minutes		
	Minimum exposure time = 20 minutes					
<b>CAUTION</b> <i>To avoid damage to the lens, do not exceed recommended exposure time.</i>						
Then:	Rinse lens <i>thoroughly</i> to remove disinfection solution. 3 cycles of 1 minute, with cool or tepid water is recommended. Dry carefully and place in a dry storage case.					
NOTE	This lens is known to be compatible with: Asepti-Wipe, Cavi-cide, Cidex, Cidex OPA, DisCide Wipe, Enviro-cide, H <sub>2</sub> O <sub>2</sub> - 3%, and Opti-Cide					
<b>CAUTION</b>	<b><i>If used on an ulcerated cornea, lens must be STERILIZED before next procedure.</i></b>					

### Sterilization

AUTOCLAVE	STERRAD	STERIS SYSTEM 1	ETO	ETO Parameters		
No	No	YES	YES	Minimum Time	Temperature	Aeration Time
		Per manufacturer instructions	See Right	1 hour	130°F (54°C)	12 hours
<b>WARNING</b> <i>Never Steam Autoclave or Boil listed lenses. Never soak in Alcohol, Acetone or Other Solvents.</i>						

For information on compatibility with alternative product care methods, contact Customer Service.



2255 116th Ave NE, Bellevue, Washington 98004-3039 USA  
 T: 425-455-5200 or 800-888-6616 F: 425-462-6669  
 E: ocular@ocular-instruments.com I: www.ocular-instruments.com

© 2002 Ocular Instruments  
 7669C3017