



## Kapetansky Water Bath

	Product Code	Size	Designed with Fredrick M. Kapetansky, M.D., Columbus, OH.	
	OKWB21 CE	21mm x 15mm x 13mm		

### Design

The saddle shape of the cup makes an ideal fit for the anterior sclera and thereby minimizes the loss of saline solution. In addition, the fluid reservoir attached to the top of the cup provides a depth of saline which is more than adequate for the ultrasonic probe to function properly.

### Technique

Since the cup has the slimmest profile possible, the user only needs to have the patient look up while inserting the cup into the inferior cul-de-sac. When the patient looks down, the upper lid will quite naturally make contact with the superior edge of the cup for a water-tight fit.

### Cleaning

Rinse: Immediately upon removal from patient's eye, thoroughly rinse in cool or tepid water.  
 Wash: Place a few drops of mild soap on a moistened cotton ball. Gently clean with a circular motion.  
 Rinse: Thoroughly rinse in cool or tepid water, and then dry carefully with a *non-linting* tissue.  
 Then: Proceed with either disinfection or sterilization instructions.

### Disinfection

SOAK IN:	<b>Glutaraldehyde</b>	OR	<b>BLEACH</b>
	2% or 3.4% aqueous solution		10% solution mixed at: 1 part bleach to 9 parts cool tepid water
	Temperature per manufacturer instructions		Recommended exposure time = 10 minutes
	Minimum exposure time = 20 minutes		
<b>Caution</b> To avoid damage to the water bath, do not exceed recommended exposure time.			

Then: Rinse *thoroughly* to remove disinfection solution.  
 3 cycles of 1 minute, with cool or tepid water is recommended.  
 Dry carefully and place in a dry storage case.

**NOTE** Material is known to be compatible with: Asepti-Wipe, Cavi-cide, Cidex, Cidex OPA, DisCide Wipe, Enviro-cide, H<sub>2</sub>O<sub>2</sub>-3%, and Opti-Cide.

### Sterilization - Autoclave

Prep: Place in sterilization case.

Process:	<b>Standard Cycle (wrapped)</b>				
	Temperature	Time	OR	Temperature	Time
	270°F (134°C)	15 minutes minimum		250°F (121°C)	30 minutes minimum

Store: Place in a biological barrier peel pouch to ensure sterility after the process.

### For Immediate Use Only

Flash autoclave (unwrapped) at a minimum of 270°F (134°C) for a minimum of 10 minutes.

<b>Sterilization - ETO</b>		
Minimum Time	Temperature	Aeration Time
1 hour	130°F (54°C)	12 hours

<b>Sterilization - Steris System 1</b>
Follow manufacturer's instructions.

For information on compatibility with alternative product care methods, contact Customer Service.



2255 116th Ave NE, Bellevue, Washington 98004-3039 USA  
 T: 425-455-5200 or 800-888-6616 F: 425-462-6669  
 E: [ocular@ocular-instruments.com](mailto:ocular@ocular-instruments.com) | [www.ocular-instruments.com](http://www.ocular-instruments.com)

© 2004 Ocular Instruments  
 8522A3017