




Ocular Pollack Iridotomy-Gonio Laser Lens

	Product Code	Contact OD	Lens Height	Lens Mag	Laser Mag	<i>Designed with:</i> Irvin P. Pollack, M.D. Baltimore, MD U.S. Pat. No. 6,698,886	
	OPIG 	15mm	21.2mm	1.5x	.65x		

Design

- The Pollack Iridotomy-Gonio Laser Lens is designed with two coated glass buttons on the anterior surface that enable performance of iridotomy and gonioscopy without changing lenses and with minimal refocusing of the slit lamp.
- Permits viewing of the anterior chamber angle.
- Designed to easily determine if the angle has been opened following iridotomy.
- The 1.5x magnification button allows lower levels of energy to be employed during the procedure.

Technique

- Apply the lens, using gonioscopy solution, in the usual manner.
- Aim the laser beam through the large plano-convex button.
- Focus sharply on the desired iris site.
- The angle can be viewed best by slight repositioning of the joystick so that the light beam falls in the angle.

Cleaning

Rinse: Immediately upon removal from patient's eye, thoroughly rinse in cool or tepid water.

Wash: Place a few drops of mild soap on a moistened cotton ball. Gently clean with a circular motion.

Rinse: Thoroughly rinse in cool or tepid water, then dry carefully with a *non-linting* tissue.

Then: Proceed with either disinfection or sterilization instructions.

Disinfecting

Soak In:	GLUTARALDEHYDE	OR	BLEACH
	2% or 3.4% aqueous solution		10% solution mixed at: 1 part bleach to 9 parts cool tepid water
	Temperature per manufacturer instructions		Recommended exposure time = 10 minutes
	Minimum exposure time = 20 minutes		
Caution <i>To avoid damage to the lens, do not exceed recommended exposure time.</i>			
Then:	Rinse lens <i>thoroughly</i> to remove disinfection solution. 3 cycles of 1 minute, with cool or tepid water is recommended. Dry carefully and place in a dry storage case.		
NOTE	This lens is known to be compatible with: Asepti-Wipe, Cavi-cide, Cidex, Cidex OPA, DisCide Wipe, Enviro-cide, and Opti-Cide.		
Caution	<i>If used on an ulcerated cornea, lens must be STERILIZED before next procedure.</i>		

Sterilizing						
AUTOCLAVE	STERRAD	STERIS SYSTEM 1	ETO	ETO Parameters		
No	No	YES	YES	Minimum Time	Temperature	Aeration Time
		Per manufacturer instructions	See Right	1 hour	130°F (54°C)	12 hours
WARNING	<i>Never Steam Autoclave or Boil listed lenses. Never soak in Alcohol, Acetone or Other Solvents.</i>					

For information on compatibility with alternative product care methods, contact Customer Service.



2255 116th Ave NE, Bellevue, Washington 98004-3039 USA
 T: 425-455-5200 or 800-888-6616 F: 425-462-6669
 E: ocular@ocular-instruments.com I: www.ocular-instruments.com

© 2001 Ocular Instruments
 7145E3017