



## Ocular Ritch Panoramic Surgical Gonioprism

	<b>Product Code</b>	<b>Image Mag</b>	<b>Contact OD</b>	<b>Handle Length</b>	<b>Static Field of View</b>	Designed with: Robert Ritch, M.D., New York, NY	
	ORPSG CE	.73x	10.8mm	88.1mm	160°		

### Lens Design

- The Ritch Panoramic Surgical Gonioprism is a glass lens designed for easy manipulation for goniotomy and direct viewing gonioscopy for adults and children.
- 180° goniotomy can be done.
- It provides a clear view of the angle during goniotomy surgical procedures.
- Lens shape creates exposed portion of cornea allowing use of instruments, incisions and corneal retraction sutures.

### Cleaning

Rinse: Immediately upon removal from patient's eye, thoroughly rinse in cool or tepid water.

Wash: Place a few drops of mild soap on a moistened cotton ball. Gently clean with a circular motion.

Rinse: Thoroughly rinse in cool or tepid water, then dry carefully with a *non-linting* tissue.

Then: Proceed with either disinfection or sterilization instructions.

### Disinfection

Soak In:	<b>GLUTARALDEHYDE</b>	<b>OR</b>	<b>BLEACH</b>
	2% or 3.4% aqueous solution		10% solution mixed at: 1 part bleach to 9 parts cool tepid water
	Temperature per manufacturer instructions		Recommended exposure time = 10 minutes
	Minimum exposure time = 20 minutes		
<b>Caution</b> <i>To avoid damage to the lens, do not exceed recommended exposure time.</i>			
Then:	Rinse lens <i>thoroughly</i> to remove disinfection solution. 3 cycles of 1 minute, with cool or tepid water is recommended.  Dry carefully and place in a dry storage case.		
NOTE	These lenses are known to be compatible with: Asepti-Wipe, Cavi-cide, Cidex, Cidex OPA, DisCide Wipe, Enviro-cide, H <sub>2</sub> O <sub>2</sub> -3%, and Opti-Cide.		

**Sterilization**

AUTOCLAVE	STERRAD	STERIS SYSTEM 1	ETO
YES	No	YES	YES

**Sterilization - Autoclave**

Prep:	Place lenses, occluder and forceps in sterilization case.				
Process:	<b>Standard Cycle (wrapped)</b>				
	Temperature	Time	<b>OR</b>	Temperature	Time
	270°F (134°C)	15 minutes minimum		250°F (121°C)	30 minutes minimum
	<b>Caution</b>	<i>Use only distilled water in the steam sterilizer.</i> If not distilled, mineral deposits from hard water (steam) will leave a cloudy film on the lens. The deposit can only be removed by regrinding and re-polishing the lens and repair costs approximate that of a new lens.			
Store:	Place in a biological barrier peel pouch to ensure sterility after the process.				
<b>For Immediate use only</b>	Flash autoclave (unwrapped) at a minimum of 270°F (134°C) for a minimum of 10 minutes.				

**Sterilization – Steris System 1**

**Follow manufacturer's instructions.**

**Sterilization - ETO**

Minimum Time	Temperature	Aeration Time
1 hour	130°F (54°C)	12 hours

**Sterilization for Autoclavable Lens Cleaning Cloth (OLCCA)**

Flash autoclave (unwrapped) at a minimum of 270°F (134°C) for a minimum of 10 minutes.

For information on compatibility with alternative product care methods, contact Customer Service.

