




Ocular Staurenghi 230 SLO Retina Lens

	Product Code 	Image Mag 	Static FOV 	Contact OD 	Lens Height 	<i>Designed with: Giovanni Staurenghi, M.D. Milan, Italy</i> <i>Patent Pending</i>	
	OSR230	.23x*	150°	19mm	69.5mm		

*relative to standard SLO image

Design

- The Staurenghi 230 SLO Retina Lens is intended for use in conjunction with a confocal scanning laser ophthalmoscope (SLO) to visualize structures of the retina and ocular fundus.
- It is optimized for use in obtaining high-resolution wide field fluorescein and indocyanine green angiography images.
- It is also effective in obtaining autofluorescence images and fundus reflectance images with green and infrared light.
- It is beneficial for diagnosis of diabetic retinopathy, peripheral retinal disorders such as heredity chorioretinal disorders, inflammatory diseases, and to document retinoschisis and retinal detachment.

Cleaning

Rinse: Immediately upon removal from patient's eye, thoroughly rinse in cool or tepid water.
 Wash: Place a few drops of mild soap on a moistened cotton ball. Gently clean with a circular motion.
 Rinse: Thoroughly rinse in cool or tepid water, then dry carefully with a *non-linting* tissue.
 Then: Proceed with either disinfection or sterilization instructions.

Disinfecting

Soak In:	GLUTARALDEHYDE		OR	BLEACH		
	2% or 3.4% aqueous solution			10% solution mixed at: 1 part bleach to 9 parts cool tepid water		
	Temperature per manufacturer instructions			Recommended exposure time = 10 minutes		
	Minimum exposure time = 20 minutes					
Caution <i>To avoid damage to the lens, do not exceed recommended exposure time.</i>						
Then:	Rinse lens <i>thoroughly</i> to remove disinfection solution. 3 cycles of 1 minute, with cool or tepid water is recommended. Dry carefully and place in a dry storage case.					

NOTE This lens is known to be compatible with: Asepti-Wipe, Cavi-cide, Cidex, Cidex OPA, DisCide Wipe, Enviro-cide, H₂O₂ - 3%, and Opti-Cide.

Caution *If used on an ulcerated cornea, lens must be **STERILIZED** before next procedure.*

Sterilization

AUTOCLAVE	STERRAD	STERIS SYSTEM 1	ETO	ETO Parameters		
No	No	YES	YES	Minimum Time	Temperature	Aeration Time
		Per manufacturer instructions	See Right	1 hour	130°F (54°C)	12 hours

WARNING *Never Steam Autoclave or Boil listed lenses.
Never soak in Alcohol, Acetone or Other Solvents.*

For information on compatibility with alternative product care methods, contact Customer Service.

