




Ocular Woldoff NA High Magnification Vitrectomy Lens

	Product Code	Image Mag	Diopter	Static FOV	Dynamic FOV	Lens Height	<i>Designed with: Herbert S. Woldoff, M.D. Phoenix, AZ</i>	
	OWIV-HMNA	0.90x	66	57°	100°	13.6mm		
								

Lens Design

- The Woldoff NA (non-autoclavable) High Magnification Vitrectomy Lens is a single-piece lens designed for clinical situations where autoclaving is either not available or not required for quick turnaround.
- It is ideal for wide angle viewing of the posterior pole.
- Its wide field provides stereopsis well beyond the area seen by a conventional flat lens.
- The high magnification and resolution create very precise depth perception.
- It provides an excellent image for delicate work around the macula such as macular hole surgery or peeling of epiretinal membranes from the macula.
- It also is the lens of choice for videotaping macular procedures.

Technique

- The lens is held on the eye by suturing one of the Landers Lens Rings to the sclera.
- After a suitable wetting agent is placed on the cornea, the lens is placed on the cornea.
- Many surgeons do not use an inverted image contact lens until the anterior third of the vitreous has been removed and a deeper image of the vitreous cannot be obtained with normal microscope observation.
- Turn off the coaxial and oblique illumination of the microscope, since this may lead to reflections from the contact lens surfaces. Check the positions of instruments repeatedly before and during the operation, as it is very difficult to recognize the patient's crystalline lens through a contact lens.
- With the microscope set at the lowest magnification and the microscope head travel at the lowest position, move the microscope to obtain a focus on the cornea. Once the lens has been placed on the eye, focus the microscope using the focusing adjustment mechanism away from the patient to focus on the image.
- It is recommended to work at the magnification where the fundus image just fills the microscope field of view.
- Be sure the lens is seated well on the cornea. If the assistant has a poor image and you find the image good (or vice versa), it is possible only one observation beam path of the microscope is receiving and transmitting a good image. Slightly shifting the lens will correct the problem.
- Keep endo-illumination as far as possible from the retina and increase illumination at its tip. This utilizes the wide-angle effect of the lens to its fullest. Light intensity at the retina will be somewhat reduced due to the distance from the retina.

Cleaning

Rinse: Immediately upon removal from patient's eye, thoroughly rinse in cool or tepid water.

Wash: Place a few drops of mild soap on a moistened cotton ball. Gently clean with a circular motion.

Rinse: Thoroughly rinse in cool or tepid water, then dry carefully with a *non-linting* tissue.

Then: Proceed with either disinfection or sterilization instructions.

Disinfection

Soak In:	GLUTARALDEHYDE		OR	BLEACH		
	2% or 3.4% aqueous solution			10% solution mixed at: 1 part bleach to 9 parts cool tepid water		
	Temperature per manufacturer instructions			Recommended exposure time = 10 minutes		
	Minimum exposure time = 20 minutes					
CAUTION		<i>To avoid damage to the lens, do not exceed recommended exposure time.</i>				
Then:	Rinse lens <i>thoroughly</i> to remove disinfection solution. 3 cycles of 1 minute, with cool or tepid water is recommended. Dry carefully and place in a dry storage case.					
NOTE	This lens is known to be compatible with: Asepti-Wipe, Cavi-cide, Cidex, Cidex OPA, DisCide Wipe, Enviro-cide, H ₂ O ₂ -3%, and Opti-Cide.					

Sterilization

AUTOCLAVE	STERRAD	STERIS SYSTEM 1	ETO	ETO Parameters		
No	No	YES	YES	Minimum Time	Temperature	Aeration Time
		Per manufacturer instructions	See Right	1 hour	130°F (54°C)	12 hours
WARNING	<i>Never Steam Autoclave or Boil this lens. Never soak in Alcohol, Acetone or other solvents.</i>					

Sterilization for Autoclavable Lens Cleaning Cloth (OLCCA)

Flash autoclave (unwrapped) at a minimum of 270°F (134°C) for a minimum of 10 minutes.

For information on compatibility with alternative product care methods, contact Customer Service.

